

# Modeling Liquid

Layering composites just got easier

Apply, shape,  
smooth,  
enjoy!



Since 1921  
100 years of Quality in Dental



# Simplify your composite works using the brush technique

Direct composites can be tricky: they sometimes stick to your instruments, can be difficult to shape or smooth just the way you want, and at times you end up with visible margin lines...

Using a brush and a modeling resin to apply composites is a true game changer.



Quick & easy application and shaping

Using a **brush wetted with a modeling liquid** makes composites less sticky, resulting in **easier and faster application & shaping**.



A smooth finish for a shorter finishing procedure

Applying the final layer using a wetted brush also enables to **smooth it efficiently** and to achieve a very nice surface. This **simplifies the finishing procedure** and saves valuable time.



Invisible for the best aesthetic results

Although slightly yellowish during application, GC Modeling Liquid becomes **fully transparent after light-curing** – guaranteeing the **best aesthetic results**.

## GC Modeling Liquid is the perfect wetting resin for the brush technique:



- Optimal **low viscosity**
- **Completely invisible** after light-curing
- Very **simple & time-saving** procedure
- 350 drops per bottle for **350 restorations!**
- Compatible with **all direct composites**
- Storage at **room temperature**



## Goodbye tricky composites

### Apply, shape, smooth, enjoy!



Dispense one drop of GC Modeling Liquid on a pad.



Use it quickly and protect it from the light with a cover.



Touch the drop with the brush in order to moisten it.



Remove the excesses thoroughly using a paper tissue.



Use a flat brush to create a palatal enamel shell.



A round brush can be useful to create the interproximal part.



Apply the dentin paste using the flat brush.



Model the dentin mamelons using a round brush.



Model the dentin mamelons using a round brush. Thanks to the brush, the final layer will have a smooth surface.

## Also optimal for posterior cases

### A round brush helps to build the cusps and refine the design of the fissures.



The brush can be easily bent for a better posterior access.



The brush can be used to apply the material cusp by cusp.



Optimal to recreate the fissures & refine anatomy.

## Get ready for easy composite layering with GC Modeling Liquid & Gradia brushes



10001380 **GC Modeling Liquid, Refill 6mL**

10000240 **GC GRADIA, Brush Round, n°1 x10**

10000241 **GC GRADIA, Brush Flat, n°1 x10**

10000242 **GC GRADIA, Brush Holder x2**



10003387 **GC Composite Modeling Kit**

Including GC Modeling Liquid,  
Gradia Brush Flat (x5), Gradia Brush Round (x5),  
Gradia Brush Holder (x2) & Mixing Pads.

## Discover great handling and aesthetics united with GC paste composites



Gradia Direct Anterior / Posterior / X



G-aëniel Anterior / Posterior



Essentia

### **GC EUROPE N.V.**

Head Office  
Researchpark,  
Haasrode-Leuven 1240  
Interleuvenlaan 33, B-3001 Leuven  
Tel. +32 16 74 10 00  
Fax. +32 16 40 48 32  
info.gce@gc.dental  
<https://www.gc.dental/europe>

### **GC UNITED KINGDOM Ltd.**

Coopers Court Newport Pagnell  
Buckinghamshire  
MK16 8JS  
United Kingdom  
Tel. +44 1908 218 999  
Fax. +44 1908 218 900  
info.uk@gc.dental  
<https://www.gc.dental/europe/en-GB>