Jet 2000

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Instructions for use

Description:

Jet 2000 is a phosphate-bonded, graphite-free precision investment for CrCo work. It was specially developed for the speed cast technique and is suitable for casting all types of dental alloy.

Casting models: Smear Vaseline on the inside of the overlapping section of the duplicating

sleeve to prevent cracking. Measure the required amount of investment and the correct concentration of expansion liquid according to the mixing ratio for casting models. Mix manually with a spatula before vacuum mixing for 30 - 40 sec. Use small vibration for pouring the contours. An instrument is useful for this. Switch off the vibrator and pour the rest of the model. When duplicating models with silicone, it is unnecessary to harden and soak the models. When duplicating models with gel, reduce the powder : liquid ratio to 100 g : 14 - 16 ml. The liquid concentration for the model and muffle should be reduced to 50 %.

Use of pressure:

Investing under pressure is not necessary due to Jet 2000's extremely homogenous grain distribution. If investing under pressure is preferred, the pressure should be released after 15 minutes so that it does not impede the setting expansion.

Mixing ratio:

100 g powder - 18 ml liquid	180 g powder - 32 ml liquid	400 g powder - 72 ml liquid
-----------------------------	-----------------------------	-----------------------------

Concentration of the expansion liquid:

Liquid concentration (1.000 ml liquid)		
Concentration	Concentrate (ml)	dist. water (ml)
50%	500	500
65%	650	350
70%	700	300
75%	750	250

Concentrations (Liquid / Water)		
100 g	180 g	400 g
9 ml / 9 ml	17 ml / 17 ml	36 ml / 36 ml
12 ml / 6 ml	21 ml / 11 ml	47 ml / 25 ml
13 ml / 5 ml	22 ml / 10 ml	50 ml / 22 ml
14 ml / 4 ml	24 ml / 8 ml	54 ml / 18 ml

Technical data:

Mixing liquid:	SILADENT type 100 expansion liquid
Mixing ratio, Powder : liquid:	100 g : 18 ml
Mixing under vacuum:	30 - 40 sec.
Working time:	4 - 5 min.
Working temperature:	approx. 22 °C (Room temperature)
Initial setting time:	7 - 10 min.
Setting expansion:	ca. 1.00 % (Typ 100 / 75 %)
Preheating temperature range for placing in the furnace:	850 - 1.000 °C
Place in the hot furnace::	15 min. after investing

SILADENT

Removing the duplicate model:

Carefully remove the investment model after approx. 20 min. Start waxing up once the model is cool enough to handle and the wax is no longer liable to melt.

Investing:

Measure and mix Jet 2000 as described above. For the investing use the same liquid concentration as for the model. A small reduction can optimize the adhesion between model and investing. Pour with minimum vibration until the pattern is completely covered. Switch off the vibrator, fill the ring and do not vibrate again. Moulds intended for speed casting should be placed on a very damp cloth so that the dehydrated model regains some moisture and achieves a better bond with the investment in the mould.

Preheating:

Exact 15 minutes after investing, place the mould with the sprue hole pointing downwards in a furnace preheated to the correct temperature (850 - 1.050°C). We recommend a maximum preheating temperature of 900°C. This temperature is perfectly adequate for both casting the mould and final thermal expansion. Moulds are ready to cast approx. 60 min. after the final temperature is reached. During the first 15 min. after inserting the moulds, the furnace fan and extractor should remain switched off and the furnace door should not be opened because of the risk of instant combustion. Jet 2000 can also be conventionally preheated at a rate of 5 - 7°C per min. Holding stages should be avoided.

Please note:

Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. SILADENT products are subjected to strict quality controls. We reserve the right to make technical changes.

Warning:

Investments contain quartz and cristobalite. Do not breathe in dust. Danger of lung diseases (silicosis or cancer). Use a suitable dust mask!

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