 **Manufacturer**
META SYSTEMS CO. LTD.
#1214-18, Sicox tower 12F, 484,
Dunchon-daero, Jungwon-gu,
Seongnam-si, Gyeonggi-do,
13229, Korea
Tel +82-31-731-7377
Fax +82-31-731-1888
<http://www.metasys.co.kr>

 **Meta Biomed Europe GmbH**
Wiesenstr. 35,
45473 Mülheim an der Ruhr,
Germany
Tel +49-208-30991910
Fax +49-208-30991999

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Australia Distribution Facility
E3Endo Pty Ltd
24 Undara Ave,
Buddina QLD 4575, Australia

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EQ-S Sonic Powered Endo Irrigation

INDICATIONS FOR USE

The EQ-S™ is used in endodontic treatment by application of sonic energy.

COMPONENT

- EQ-S Handpiece
- EQ-S tip
- Protective Sleeve

CONTRAINDICATIONS

- The unit should not be used on a patient with a pacemaker or ICD.
- Interference may occur when used on patients with a pacemaker or ICD.
- Ask patients and users if they have an implanted device before using this product. Explain the circumstances to them.
- Weigh the risks and benefits and contact your patient's cardiologist or appropriate qualified healthcare professional prior to performing the treatment.
- Symptoms including a raised heartbeat, irregular pulse, and dizziness may signal problems with a pacemaker or ICD. Make appropriate emergency provisions and take immediate action if patients become ill.

WARNING

- Do not submerge the handpiece in water or any disinfecting solution.
- Do not sterilize the handpiece.
- The EQ-S tip is intended for a single patient use only. Do not re-use, due to risk of cross-contamination
- Do not clean, disinfect, or sterilize the EQ-S tip.

CAUTION

- U.S. federal law restricts this device to the sale by, or on the order of a dentist.
- Inspect the tip before use. If the product is damaged, do not use the product and discard it.
- To change the direction of the tip, take the tip out from the handpiece and re-attach it with the desired direction. Otherwise, rotating the tip while it is attached in the handpiece may cause the damage of the handpiece.

ADVERSE REACTIONS

- None known
- In the event of a serious accident related to the use of this medical device, report it to the manufacturer and the

competent authority of the country where the user and/or patient resides.

INTENDED PATIENT POPULATION

- Age: From children to adults
- Weight/Nationality/Sex: Not applicable

THE NAME OF EACH PART

- EQ-S tip: Vibration tip
- POWER/MODE button:
 - Turn on: Press the POWER/MODE button briefly.
 - Select desired mode: After turning on the power, press the POWER/MODE button briefly. (Amber LED: High speed mode, Blue LED: Low speed mode)
 - Turn off: Press and hold the POWER/MODE button.
- Battery cover: Rotate the battery cover counterclockwise to disengage the battery cover and remove it from the handpiece.

BEFORE USE

- This device is not supplied with the battery inserted. When using the device for the first time, please insert "AA" alkaline battery.
- Tips are packaged individually. Remove the sealed wrapping paper before use.

INSTRUCTIONS FOR USE

- Fill the pulp chamber with NaOCl and/or EDTA solution such as MD-Cleanser.
- Select the proper EQ-S tip, attach to the handpiece, and place it into the root canal.
- Insert the EQ-S tip 5mm shorter than the measured working length. (Tip length indication: 18mm and 20mm from the end of the tip)
- To activate the device, push power/mode button.
 - High speed will activate by default upon activation so use same power/mode button to control the speed level. (high or low options available)
- Move the handpiece up and down approximately 2-3mm.
- Irrigate for 30-60 seconds then use the suction to remove debris.
- Repeat above steps as needed.

EQ-S TIP REMOVAL

- Point the tip away from patient or personnel and pull it off.
- Discard used tips.

BATTERY REPLACEMENT

- Rotate the battery cover counterclockwise to open.
- Remove and replace with a new AA alkaline battery.
- Turn the battery cover clockwise to securely close.

CLEANING

Clean all the surfaces with a cloth lightly moistened with other intermediate level disinfectant.

DISPOSAL

Please contact your nearest representative office or your local environmental office in case of disposal of this unit.

SPECIFICATIONS

- High speed mode (Amber): 9,000 cpm
- Low speed mode (Blue): 8,000 cpm
- Rating: 1.5 Vd.c. (AA Alkaline battery)
- Dimension: ~183mm (Length)
- Weight: 66g (Including battery)
- According to the degree of protection against electric shock: Type B Applied part
- According to the degree of protection against harmful ingress of water: IPX0

SYMBOLS

	Serial number
	Manufacturer
	Date of manufacture
	This side up
	Fragile, handle with care
	Keep dry
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation
	Do not re-use
	CTP mark
	Medical device
	Authorized representative in the European Community/European Union
	CE mark
	CSA mark
	Warning
	Caution
	See Instructions
	Type B applied part

	WEEE wheeled bin
	Catalogue number
	Importer

OPERATING, STORAGE & TRANSPORT

	Operating conditions	Storage and transport conditions
Temperature	10 – 40°C	-20 – 60°C
Humidity	30 – 75 %	5 – 90 %
Atmosphere	700 – 1060hPa	700 – 1060hPa

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QUALITY ASSURANCE CERTIFICATE

Meta Systems Co., Ltd. guarantees device quality by offering services free of charge for failures caused by defects in the device during the warranty period.

Product name	Sonic Powered Endo Irrigation	Model	EQ-S
Warranty	Handpiece: 1 year		
Manufacturer	• Business name: Meta Systems Co., Ltd. • Address: #1214-18, Sicox tower 12F, 484, Dunchon-daero, Jungwon-gu, Seongnam-si, Gyeonggi-do, 13229, Korea • Phone: 82-31-731-7377		
Sales place info	• Business name: • Address: • Phone:		

Buyer info

Business name		Phone	
Address			
Date of manufacture		Serial number	
Date of purchase		Purchase price/quantity	

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