



ORO CLEAN® Plus

Gentle disinfection and cleaning of all types of suction units



- ▶ **Effective disinfection**
- ▶ **Non-foaming**
- ▶ **Biodegradable**

ORO CLEAN® Plus is a highly effective concentrate for the disinfection and cleaning of dental suction units, spittoon bowls, and amalgam separators. An efficient combination of disinfecting and cleaning agents successfully inactivates germs, removes debris, and ensures the smooth running of your aspirating unit. The aldehyde- and alcohol-free formula provides extensive protection of the unit and is easily biodegradable. ORO CLEAN® Plus's non-foaming action is non-damaging to the suction pump.

Use



For disinfection of suction units, use an aqueous solution and dosage as required. To measure the correct dosage, use lid as a measuring cup (lower scale = 20 ml, upper scale = 40 ml). Rinse with water prior to disinfection. Once a week and as required, aspirate 1.5 liter of working solution per unit using the OROCLEAN® Cup. Pour 0.5 liter down the spittoon bowl. After disinfection, rinse the unit thoroughly with water. In addition, rinse the unit thoroughly with water after surgical / invasive procedures.

Scope of Application



- ▶ For the disinfection and cleaning of suction units, spittoon bowls, and amalgam separators

Precautions



R36 / 38 Irritating to eyes and skin. S2 Keep out of the reach of children. S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. S28 After contact with skin, wash immediately with plenty of water. S37 / 39 Wear suitable gloves and eye / face protection. S45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). Use disinfectants safely. Always read label and product information before use.

Delivery Forms



2 liter, 5 liter



Effective Spectrum / Exposure Times

Bactericide	2% in 30 min.
Fungicide (<i>C. albicans</i>)	2% in 30 min.
Hepatitis B / HIV (1)	2% in 60 min.
BVDV / Hepatitis C (2)	2% in 5 min.

1) HBsAg antigen test; 2) BVDV surrogate virus test

Complies with DGHM Standard Methods and EN 1040, EN 1275, AFNOR NF T 72150, NF T 72170



Dosing

	Concentrate + Water	Time
2%	40 ml + 1960 ml	60 min.

Composition in 100 g

5 g Alkylbenzyltrimethyl ammonium chloride, Antifoaming agents, Special surfactants.