




Safety visor properties



No
distortion



Anti-fog
Resistant to repeated washing



FEATURES

- > Interchangeable visor in acetate
- > Adjustable elastic band
- > Adjustable inclination angle
- > No optical distortion (Optical Class 1)
- > Anti-fog - Resistant to repeated washing (N marking)
- > Impact resistant (S marking)

CODE

MDU7015S
Safety visor

PROTECTION

Facial protection
Certified according to EN 166 point 3 for drops and splashes

COMPATIBILITY WITH
MAGNIFYING SYSTEMS & HEADLIGHTS

Usable with Galilean system	•
Usable with Prismatic system	•
Usable with headlight EOS 2.0	•
Usable with headlight EOS HP	•
Usable with headlight EOS Wireless	•

MARKING

U 1 S N CE - U 166 3 S CE
U Z87 X - U Z87 D3