

The simple solution for all cavities

Evetric® Line

Your everyday composites

NEW



The simple solution for all cavities

The Evetric® Line now comprises four products to cover all your composite needs: Evetric®, Evetric® Flow, Evetric® Bulk Fill and Evetric® Bulk Flow.

Are you looking for a flowable consistency? Do you want to be more efficient in the posterior region? The Evetric Line has everything covered for you.

The core of the line is Evetric®

Evetric – your trusted composite – is a sculptable universal composite for the anterior and posterior region.

- ✓ Lifelike shade blend due to the pronounced chameleon effect and different translucencies for enamel and dentin shades
- ✓ Composite with high flexural strength^[1] of 120mpA and high radiopacity
- ✓ Easy handling and long processing time (more than 3 minutes)
- ✓ Easy to contour as a result of low stickiness, good spreadability and creamy consistency.
- ✓ Efficient polishability for a lustrous surface

Recommendation for use

- ✓ Use universally applicable enamel shades for small defects and direct veneers with a translucency of 11.5%
- ✓ Use dentin shades showing low translucency (7%) for restoring large defects and fulfilling high esthetic expectations

[1] According to ISO 4049

Exact application with Evetric® Flow

Evetric Flow is the flowable universal composite for tooth-preserving restorative treatments in the anterior region. It can also be used as a liner in combination with Evetric.



- ✓ Lifelike shade blend due to the pronounced chameleon effect
- ✓ Composite with high flexural strength^[2] of 120mpA and high radiopacity
- ✓ Easy handling and long processing time (more than 3 minutes)
- ✓ Precise application due to thixotropic consistency
- ✓ Efficient polishability for a lustrous surface
- ✓ Proven, time-tested quality

[2] According to ISO 4049

Easy modelling with Evetric® Bulk Fill

Evetric Bulk Fill is a sculptable bulk composite designed to deliver natural-looking esthetic results in the posterior region.

- ✓ High esthetics for posterior restorations thanks to patented technologies
- ✓ Reliable curing of increments of up to 4 mm in thickness with curing time of 5 seconds
- ✓ High flexural strength, low shrinkage and low susceptibility to process-related air entrapments
- ✓ Effective combination with Evetric® BulkFlow
- ✓ Easy to distinguish on X-rays due to high radiopacity
- ✓ Long processing time (more than 4 minutes)

Fast volume build-up with Evetric® Bulk Flow

Evetric Bulk Flow is a flowable bulk composite designed to deliver natural-looking esthetic results in the posterior region.

- ✓ High esthetics for posterior restorations thanks to patented technologies
- ✓ Reliable curing of increments of up to 4 mm in thickness with curing time of 5 seconds
- ✓ High flexural strength, low shrinkage and low susceptibility to process-related air entrapments
- ✓ Effective combination with Evetric BulkFill and Evetric



Fewer treatment steps due to 4-mm increments

Compared to conventional filling techniques, cavities up to 4 mm in depth can be restored with a single increment. For cavities deeper than 4 mm, increments up to 4 mm thick can be applied. Using fewer increments reduces the need to constantly switch between the composite, modeling instrument and light-curing device, while still delivering the same high-quality treatment results as conventional layering techniques.

From conventional layering to bulk-fill technique



Highlight

Excellent for treating deep cavities thanks to self-levelling effect and great adaptation to cavity walls

Your benefits with the Evetric Line at a glance

Four composites are all you need for an efficient restorative workflow for all cavity classes^[3].

Lifelike esthetics

Unlock lifelike esthetics and effortless anatomical contouring in the anterior and posterior region with our products for direct restorative treatments.

High-quality results

High quality is our priority. The products are precisely coordinated to work seamlessly together, enabling efficient and reliable application and processing.

Efficiency and easy

The bulk fill materials are easy to use thanks to the standardized curing times and have increased efficiency due to the 4 mm layers.

The winning team for successful direct restorative treatments

Evetric® Bond and Bluephase® EasyCure complement the Evetric composites and will allow you to focus on what really matters: your patients.



Evetric Bond is the state-of-the-art self-etch adhesive for convenient bonding. A single-component adhesive with high bond strength and 10 seconds curing time



Bluephase® EasyCure is the easy-to-use, everyday curing light for effortless polymerization. A product that stands out for its quality and simplicity. With the press of a single button, you can activate the High-Power program for fast, reliable and efficient curing in just 10 seconds – every time.

Esthetics, quality and simplicity with well-coordinated products

	Area of application	Our recommendation for cavity classes*	Consistency	Shade	Curing time with Bluephase EasyCure	Curing time with Bluephase PowerCure
Evetric® Bond	Adhesive / Bonding agent	–	–	–	10 seconds at 1200 mW/cm²	10 seconds at 1200 mW/cm²
Evetric®	Filling material, max. layer thickness 2 mm	3, 4 and 5	Sculptable	A1 / A2 / A3 / A3.5 / B2 / B3 / C2 / A3.5 Dentin / B2 Dentin / T	10 seconds at 1200 mW/cm²	10 seconds at 1200 mW/cm²
Evetric® Flow	Filling material, max. layer thickness 2 mm	3, 4 and 5	Flowable	A1 / A2 / A3 / A3.5	10 seconds at 1200 mW/cm²	10 seconds at 1200 mW/cm²
Evetric® Bulk Flow	Filling material, max. layer thickness 4 mm	1 and 2	Flowable	IVA	10 seconds at 1200 mW/cm²	5 seconds at 2100 mW/cm²
Evetric® Bulk Fill	Filling material, max. layer thickness 4 mm	1 and 2	Sculptable	IVA	10 seconds at 1200 mW/cm²	5 seconds at 2100 mW/cm²

[3] Cavity classes I–V according to G. V. Black

* Further areas of application and information can be found in the IFU.

Direct restoratives for efficient esthetics

Our well-coordinated products for direct restorative procedures enable you to efficiently create high-quality, esthetic restorations.



Isolate

Relaxed and efficient treatments with OptraGate 2



Bond

Powerful bonding with Evetric® Bond



Restore

Evetric® Line – the simple solution for all cavities



Cure

High curing performance delivered by Bluephase EasyCure



Polish

Lustrous surfaces produced with OptraGloss



Protect

Effective prevention with Fluor Protector S