

THE PREMIUM PROPHYLAXIS SYSTEM

# LUNOS®

MAKING SMILES BRIGHTER



## LUNOS® HAS FOUR REASONS FOR BRIGHTER SMILES

- 1.** Systematic prophylaxis. Fast, effective, and convenient in application. For all prophylaxis specialists, dentists, hygienists and patients
- 2.** Premium prophylaxis solution with an extensive product portfolio ranging from handpieces to powder and pastes
- 3.** Meets the highest medical standards in terms of product performance and ingredients
- 4.** Low-pain treatment and maximum patient comfort



 HENRY SCHEIN®  
1300 65 88 22 | henryschein.com.au

 DÜRR  
DENTAL  
THE BEST, BY DESIGN

## EASY TO HANDLE, VERSATILE, WITH RADIANCE.

The powder jet handpiece MyLunos® offers simple, powerful and ergonomic assistance when you need to remove discolouration, deposits and biofilms. Using one single handpiece with two different nozzles, it is ready for either supragingival or subgingival application. The detachable powder chamber is quick and easy to change, leaving one free to always use exactly the right powder for the purpose. MyLunos® offers great flexibility and efficiency.

### HANDPIECE AND INDIVIDUAL PARTS CAN BE REPROCESSED FULLY AUTOMATICALLY

All individual elements of the handpiece are made of temperature-resistant plastic. And can therefore be thermo-disinfected and autoclaved.




### EXCHANGEABLE CHAMBER PRINCIPLE

Easy to change, even during treatment. No work interruptions or loss of power due to a low fill level. Different powders can also be used to suit requirements.



### 360 DEGREE ROTATABLE TIP

The tip can be rotated by 360 degrees enabling particularly ergonomic work.



### FOR SUBGINGIVAL TREATMENT

The Perio tip is delivered sterile packaged and can be safely disposed of after application.

(Working length approx. 11.7mm)

## POWERFUL YET GENTLE: LUNOS® POWDERS

- Innovative abrasive agent trehalose
- Best water solubility
- Ergonomic bottle with a single-hand flip-spout cap



### LUNOS® PROPHY POWDER GENTLE CLEAN

For supragingival, fast and effective cleaning and removal of extrinsic discolouration.

Flavours: Neutral, Orange and Spearmint. 180g bottle

### LUNOS® PROPHY POWDER PERIO COMBI

For subgingival and supragingival cleaning. Ideal for recall. Removes biofilm quickly, effectively, thoroughly and yet gently.

Flavour: Neutral. 180g bottle

## POWERFUL YET GENTLE: LUNOS® PASTES



### LUNOS® PROPHY PASTE TWO IN ONE

The innovative two-in-one solution: removes discolouration and polishes in a single step. Saves time, yet is thorough and gentle on tooth substance.

The abrasiveness of the paste decreases gradually thanks to the self-reducing polishing particles

- Can be used for all treatment needs
- Effective cleaning and polishing thanks to an innovative system of abrasive and polishing agents. Hydroxyapatite supports remineralisation
- Particularly gentle on tooth substance, RDA 28
- Mint flavour additionally contains sodium fluoride to support remineralisation

Flavours: Neutral and Mint. 100g tube



### LUNOS® PROPHY PASTE SUPER SOFT

Lasting, gentle fine polish, even on particularly sensitive surfaces.

- Universal use, even on implants, orthodontic appliances and dental restorations
- Creates smooth, high-gloss surfaces, and new layers of plaque or calculus take longer to form
- Hydroxyapatite cleans and polishes extremely gently and supports remineralisation
- RDA < 5

Flavours: Neutral and Orange 50g tube



### LUNOS® PROPHYLAXIS RING

For optimised handling through simple application of the propy paste.

- Flexible ring size
- Re-usable



### LUNOS® PROPHYLAXIS CUSHION

Provides safe and comfortable head support. Supports head and neck thanks to viscoelastic foam material.

- Adapts to each patient
- Washable and disinfectable

Dürr Dental SE  
PO Box 2067  
Woonona East NSW 2517  
Australia  
email: [info@duerrdental.com](mailto:info@duerrdental.com)  
[www.duerrdental.com/en/AU](http://www.duerrdental.com/en/AU)



Scan QR Code for  
product information

**HENRY SCHEIN®**  
1300 65 88 22  
[henryschein.com.au](http://henryschein.com.au)

**DÜRR DENTAL**  
THE BEST, BY DESIGN