



*... and everything to
make you successful.*

Dealer information

Charisma Bulk Flow ONE



Product Description



CHARISMA BULK FLOW ONE - ONE SHADING MEETS BULK FILLING.

Charisma Bulk Flow ONE simplifies and streamlines the restorative process

As a further universal shade solution for efficient restorations in the posterior region, it represents the first flowable **one-shade bulk-filling composite** in the dental market that adapts to any tooth shade from A1 to D4.



- **Always grab the right shade:** With the incredible shade-matching properties of Charisma Bulk Flow ONE the restoration becomes invisible by blending into the surrounding dentition.



- **Efficient handling:** Cavities are filled quickly and easily with increments of up to 4 mm without the need for an additional capping layer.



- **Long-lasting restorations:** Good margin adaptation = tight margins, low shrinkage stress, high flexural strength and high wear resistance ensure durable restorations that withstand the high mastication load in the posterior region.

What makes **Charisma Bulk Flow ONE** unique?

It is the first flowable bulk fill composite that can be used without additional capping.



Source: Prof. Dr. Claus-Peter Ernst, Mainz Germany

Light-curing times and maximum layer thickness **Charisma Bulk Flow ONE:**

- Curing depth 4 mm
- Curing time 20 sec

Indications for **Charisma Bulk Flow ONE:**

- Class I, II, III and V direct restorations
- First layer lining for class I and II cavities
- Repair of direct and indirect restorations in combination with a suitable adhesive
- Splinting of loosened teeth resulting from trauma or periodontal associated events
- Extended fissure sealing.

Selling units & pricing: Refills and Value-Packs

Charisma Bulk Flow ONE - Refills			Recommended retail price
66094725	CHARISMA BULK FLOW ONE SYR 1 x 2G	1 x 2g syringe 5 x cannulas 1 x IFU	29,90 € Purchase price 20.03 €

Charisma Bulk Flow ONE – Value Packs			Recommended retail price
66095336	CHARISMA BULK FLOW ONE VALUE-KIT SYRINGE	10 x 2g syringes 10 x 5 cannulas 1 x pictorial card, 1 x IFU	249 € Purchase price 166.83 €

-17%

