



elements™ Connect & Apex Connect

Cordless Endodontic Motor and Electronic Apex Locator

Seamless Connection. Greater Confidence.

Add confidence to your shaping, with intuitive elements Connect and Apex Connect. Designed to work seamlessly together, interconnected system supports patient safety while simplifying your procedures.

elements™ Connect and Apex Connect

An Interconnected System



Dual Light Indicators

With intuitive light indicators in both motor handpiece and the Apex file clip, you'll have greater visibility for more confident shaping.



Patient Safety

The intuitive, easy-to-see smart light and audible indicators, as well as a large, bright viewing display alert you when you're reaching the apex. May be used independently or as a connected system.



Portable Convenience

Easily move between operatories with this interconnected system; lightweight, ergonomic and cordless design.



Versatile Motions

With Adaptive, Reciprocation and Rotary motions in one motor, elements Connect allows offices to maintain their existing techniques and protocols.


Preset

9 Kerr preset file settings, 1 reciprocating setting and 4 customizable settings.

Continuous Use



Two removable batteries are included to ensure one is always charged so you can complete your procedures safely without interruption.



Built-in Accuracy

Ensure measurement accuracy with Apex Connect's automatic self-calibration feature that activates every time the device is turned on.

Shape with Confidence

View real-time working length information during procedures with a smart light indicator that changes color, as well as audible indicators, as you reach the apex.

Ordering Information

elements Connect

- 815-1800 elements Connect Cordless Motor Kit
- 815-1806 elements Connect AC/DC Adapter
- 815-1808 elements Connect Contra Angle
- 815-1809 elements Connect Charging Base
- 815-1810 elements Connect Battery Pack

Apex Connect

- 815-1802 Apex Connect Kit - US
- 815-1807 Apex Connect AC/DC Adapter
- 815-1812 Apex Connect Probe Cord
- 815-1813 Apex Connect Lip Hook
- 815-1814 Apex Connect File Hold Gender
- 815-1815 Apex Connect File Holder A
- 815-1816 Apex Connect File Holder B
- 815-1817 Apex Connect Replacement Battery

Power Cord

- 815-1804 elements/Apex Connect Power Cord



**To Learn More About Kerr
Endodontic Solutions
Contact Your Kerr Rep**



**For More Information About elements Connect
and Apex Connect Scan the QR Code**